



Absorbable Surgical Sutures Catgut Sutures-Plain

Stargut, a natural absorbable suture is obtained by longitudinally slitting the intestinal submucosa of sheep and goat, twisting the ribbons and joining them in wet condition giving monofilament finish.

Monofilament

For smooth tissue passage. Smooth tie down and knot security

Tissue Reaction

Moderate. Plain catgut enjoys lower tissue reaction as compared to chromicised

Absorbable

By enzymatic digestion. Complete by 90 days.

Tensile Strength

Strength retention for at least 21 days.

Colour

Straw

Range

Gauge sizes 5/0 to 2 (U.S.P) Supplied as needled sutures and ligatures.

Sterilisation

Ethylene Oxide.

Caution

Avoid using where long term tissue approximation is needed.
Absorbtion is faster in infected tissues.

Frequent Applications

General Closure, Obstetrics And Gynaecological Surgeries Ophtalmology, Urology, Orthopaedics, Bowel Anastomosis Ligation







STERILISED SURGICAL CATGUT SUTURE: PLAIN

EYELESS NEEDLED SUTURES

NEEDLE	Needle Profile	SUTURE	Suture Size in USP				
LENGTH	& Description	LENGTH	5-0 (1.5)	4-0 (2)	3-0 (3)	2-0 (3.5)	
12mm	3/8 Circle Reverse Cutting	76 cm	FN 4269				
16 mm	3/8 Circle Cutting	76 cm		FN 4270	FN 4271		
22mm	3/8 Circle Cutting	76 cm			FN 5010		
25mm ⑤	½ Circle Rounnd Bodied	76 cm			F M 238		
26 mm	3/8 Circle Cutting	76 cm			FN 5020		
30 mm (i)	3/8 Circle Rounnd Bodied	76 cm				F ⊠ 033	
50mm	Straight Rounnd Bodied	76 cm			F ⊠ 007	F15100 8	
60mm	Straight Cutting	76 cm			FN 4288		

STANDARD AND SHORT LENGTHS

SUTURE LENGTH	Suture Size in USP								
	5-0 (1.5)	4-0 (2)	3-0(3)	2-0 (3.5)	0 (4)	1 (5)	2 (6)	3 (7)	4 (8)
1 X 152 cm	FS 2000	FS 2001	FS 2002	FS 2003	FS 2004	FS 2005	FS 2006	FS 2007	FS 2008
2 X 76 cm	FS 2240	FS 2241	FS 2242	FS 2243	FS 2244	FS 2245	FS 2246	FS 2247	FS 2248